

Ningbo Dasky Life Science Co., Ltd



Address: No.6 Donggong One Road, Dongguang Industrial park, Jiangshan, Ningbo

Tel: +86 574 88455466 Fax: +86 574 88459578

Email: sales@daskymed.com

Website: www.daskymed.com

STATEMENT

To Whom It May Concern:

We Ningbo Dasky Life Science Co., Ltd, here by certify that the above produce manufactured by us conforms to the following points:

MEDIUM COMPONENTS

Hanks' Balanced Salts
Sucrose
Penicillin
Gentamicin
Streptomycin sulfate
Amphotericin B
Nunessential amino acid
Phenol Red
Ultrapure water



VSM01 medium ingredients of hanks balanced salt as following points:

Components	CAS #.	% by wt.
Nacl	7647-14-5	0.8%
MgSO4•7H2O	10034-99-8	0.01%
Kcl	7447-40-7	0.04%
Mgcl2•6H2O	7786-30-3	0.01%
Cacl2	10043-52-4	0.014%
Na2HPO4•12H2O	10039-32-4	0.0152%

KH ₂ PO ₄	7778-77-0	0.006%
Glucose	50-99-7	0.1%
Ultrapure water		99.0048%

The Nonesential amino acid in VSM01 medium components as following points:

Component	Molecular weight	Concentration	Molarity (mM)
Glycine	75	750	10
L-Alanine	89	890	10
L-Asparagine	132	1320	10
L-Aspartic acid	133	1330	10
L-Glutamic Acid	147	1470	10
L-Proline	115	1150	10
L-Serine	105	1050	10

COMPANY NAME: Ningbo Dasky Life Science Co., Ltd

POSITION: Quality Control

SIGNATORY: WANG XIANWEI

Date: 2021.03.19

